

ActiV.A.C.[®] Therapy System

Designed Specifically for the Ambulatory Patient

Introducing the ActiV.A.C.[®] Therapy System, a lightweight, portable system designed specifically for the ambulatory patient. It combines enhanced technology with easy-to-use new features to facilitate the mobility of patients who may benefit from V.A.C.[®] Therapy.

LIGHTWEIGHT, PORTABLE THERAPY UNIT

- 2.4 lbs. – 34% lighter*
- 14-hr average battery life – 17% longer*
- 36% smaller*
- Discreet – new carrying case hides tubing
- Informative – enhanced alarm identification*
- Efficient – easy, quick release 300 mL canister
- USB port for convenient data download

SENSA.T.R.A.C.[®] TECHNOLOGY

- Helps monitor and maintain target pressure
- Alarm differentiation for easier troubleshooting
- Seal Check™ to identify and resolve leaks

SENSA.T.R.A.C.[®] PAD

- Designed with patient comfort in mind
- More flexible pad material for easier application over body contours
- Designed with enhanced fluid dynamics to help reduce tubing blocks and associated alarms
- Low profile design is discreet under clothing

SEAL CHECK™

- Helps clinicians and patients identify and troubleshoot negative pressure leaks
- Provides real-time feedback with audible and visible cues for leak location
- May reduce unscheduled clinician visits

THERAPY HISTORY REPORT

- Allows clinicians to document compliance with V.A.C.[®] Therapy

*Compared to the V.A.C. Freedom[®] Therapy System

ActiV.A.C.® Therapy System Specifications

Weight:2.4 lbs (1.08 kg) with an empty canister
 Dimensions:7.6" W x 6" H x 2.5" D
 (19.3 x 15.2 x 6.4 cm)
 Battery Type:Lithium ion rechargeable battery
 Battery Life:14 hour average
 Recharging Time:6 hours to reach full charge from a
 completely discharged state
 Canister Size:300 mL
 Therapy Delivery Modes:Continuous or Intermittent
 Intensity Settings:Low, Medium, High
 Therapy Pressure Options:Negative Pressure: -25 mmHg through
 -200 mmHg (-3.3 kPa to -26.6 kPa)
 Carrying Case:The black nylon design allows for full
 touch screen navigation, battery
 charging, canister replacement and
 tubing storage without removing the
 ActiV.A.C.® Therapy Unit from the case.
 The case has an adjustable, padded
 strap and a belt loop that allows the
 unit to be worn in multiple orientations.
 A hook is included for hanging.

Electrical Data

External Power Supply Input:100-240 VAC .72 A @ 115 VAC 47-63 Hz
 External Power Supply Output:12 V, 3.3 A
 Patient & Enclosure Leakage Current: <100 Microamps

IEC Classification

- Type B, Applied Part Class II IPX0
- Equipment not suitable for use in presence of flammable anesthetic mixture with air, oxygen or nitrous oxide

Storage Conditions

Temperature Range:-4°F (-20°C) to 140°F (60°C)
 Relative Humidity Range:0-95% non-condensing

Operating Conditions

Temperature Range:41°F (5°C) to 104°F (40°C)
 Relative Humidity Range:0-95% non-condensing
 Altitude Range:0 to 14,000 feet (0 to 4267 m)
 Optimum Performance:0 to 8,000 feet (0 to 2438 m)

Ordering Information

Part Number	Description	Case Quantity
V.A.C.® Dressings with SensaT.R.A.C.® Pad		
M8275051/10	V.A.C.® GranuFoam® Small Dressing	10
M8275051/5	V.A.C.® GranuFoam® Small Dressing	5
M8275052/10	V.A.C.® GranuFoam® Medium Dressing	10
M8275052/5	V.A.C.® GranuFoam® Medium Dressing	5
M8275053/10	V.A.C.® GranuFoam® Large Dressing	10
M8275053/5	V.A.C.® GranuFoam® Large Dressing	5
M8275098/10	V.A.C. GranuFoam Silver® Small Dressing	10
M8275098/5	V.A.C. GranuFoam Silver® Small Dressing	5
M8275096/10	V.A.C. GranuFoam Silver® Medium Dressing	10
M8275096/5	V.A.C. GranuFoam Silver® Medium Dressing	5
M8275099/10	V.A.C. GranuFoam Silver® Large Dressing	10
M8275099/5	V.A.C. GranuFoam Silver® Large Dressing	5
M8275064/5	V.A.C.® GranuFoam® Hand Dressing	5
M8275065/5	V.A.C.® GranuFoam® X-Large Dressing	5
M6275033/10	V.A.C.® WhiteFoam Small Dressing Foam Only 10 x 7.5 cm	10
M6275034/10	V.A.C.® WhiteFoam Large Dressing Foam Only 10 x 15 cm	10
M8275067/10	V.A.C.® WhiteFoam Large Dressing 10 x 15 cm	10
M8275067/5	V.A.C.® WhiteFoam Large Dressing 10 x 15 cm	5
M8275068/10	V.A.C.® WhiteFoam Small Dressing 10 x 7.5 cm	10
M8275068/5	V.A.C.® WhiteFoam Small Dressing 10 x 7.5 cm	5

Part Number	Description	Case Quantity
V.A.C.® Dressings with SensaT.R.A.C.® Pad		
M8275074/5	V.A.C.® GranuFoam® Heel Dressing	5
M8275075/10	V.A.C.® GranuFoam® Round Dressing	10
M8275075/5	V.A.C.® GranuFoam® Round Dressing	5
M8275081/10	V.A.C.® GranuFoam® Thin Dressing	10
M8275081/5	V.A.C.® GranuFoam® Thin Dressing	5
M8275080/5	V.A.C.® Abdominal Dressing System with Non-Adherent Layer	5

Accessories

M8275057/10	SensaT.R.A.C.® Pad Only	10
M6275009/10	V.A.C.® Drape	10
M6275097/5	V.A.C.® Extra Large Drape	5
M6275066/10	T.R.A.C.® "Y" Connector	10
M6275066/5	T.R.A.C.® "Y" Connector	5
M6275069/10	T.R.A.C.® Tubing Cap	10
M6275069/5	T.R.A.C.® Tubing Cap	5

ActiV.A.C.® Canisters

M8275058/10	300 mL ActiV.A.C.® Canister w/Gel (Sterile)	10
M8275058/5	300 mL ActiV.A.C.® Canister w/Gel (Sterile)	5



ActiV.A.C.® Therapy System
Part # 340000

To learn more about the new
 ActiV.A.C.® Therapy System, visit
www.ActiVAC.com or call 1-800-275-4524.

